

NOVEMBER 2013

NEW SALES ARGUMENT FOR HELIPROBE!

Comparison study between Heliprobe[®] System (HeliCap[™]) and Py-test (PYCap)

A study performed by Prof Tom Borody et al, Sydney Australia, was presented at the DDW (Digestive Diseases Week) in Orlando, 18th-22nd May. Patients were tested with the two UBT-methods and endoscopy (histology and RUT).

The conclusion from the study was that HeliCap[™] is superior to PYCap, regarding all parameters analysed, such as sensitivity, specificity, NPV and PPV. HeliCap[™] diagnosed the infection with an accuracy of 95%, considerably greater than the 83,3 % seen with PYCap.

Please use this study, enclosed with this Newsletter, in discussion with customers. The Heliprobe[®] System is superior to PY-test and also much more easy to use.

The Heliprobe[®] System only requires a capsule to be swallowed and then the patient exhales into a BreathCard[™] which is analysed by the Heliprobe[®] instrument. When performing the PY-test you need scintillation liquids and a scintillation counter, which is not very user-friendly.

Best regards
KIBION